


Lecture 1

Introduction to the Concept of Virtual Instrument (VI) and Graphical Programming Language

吳文中

What Is LabVIEW?

— A graphical programming environment used to develop sophisticated measurement, test, and control systems.



LabVIEW:


- Interfaces with wide variety of hardware
- Scales across different targets and OSs
- Provides built-in analysis libraries

Graphical Language


- **Why NI LabVIEW?**

- **Graphical Programming for Test, Measurement, and Control**
- Rapid application development with Express VIs and easy-to-use graphical environment
- Interactive measurement assistants and powerful redesigned DAQ interface for connecting to all types of I/O
- Expanded targeting options from Real-Time to FPGA



- **Competitors? HP(Agilent) VEE**



A History of LabVIEW




- **Modern LabVIEW ~**
- **August 2005** • LabVIEW 8 Project-based developments, distributed intelligent
- **May 2003** • LabVIEW 7 Express VIs, I/O Assistants, FPGA/PDA targets
- **January 2002** • LabVIEW 6.1 Enhanced networking capabilities, analysis
- **August 2000** • LabVIEW 6/Internet-ready measurement intelligence
- **March 1998** • LabVIEW 5.0 ActiveX, Multithreading
- **February 1996** • LabVIEW 4.0 Added professional tools, improved debugging
- **August 1993** • LabVIEW 3.0 Multiplatform version of LabVIEW
- **September 1992** • LabVIEW for Windows
- **January 1990** • LabVIEW 2.0 for Macintosh
- **October 1986** • LabVIEW 1.0 for Macintosh
- **April 1983** • LabVIEW project begins



Acquire, Analyze, and Present

Nearly all test, measurement, and control applications can be divided into 3 main components: the ability to *acquire*, *analyze*, and *present* data. LabVIEW is the easiest, most powerful tool for acquiring, analyzing, and presenting real-world data.

Acquire with LabVIEW




LabVIEW is tightly integrated with all NI hardware, in addition to connecting to thousands of I/O devices from hundreds of different vendors.

LabVIEW can acquire data using the following devices and more:

- GPIB, Serial, Ethernet, VXI, PXI Instruments
- Data Acquisition (DAQ)
- PCI eXtensions for Instrumentation (PXI)
- Image Acquisition (IMAQ)
- Motion Control
- Real-Time (RT) PXI
- PLC (through OPC Server)
- PDA
- Modular Instruments

Analyze with LabVIEW



Powerful measurement analysis is built in to the LabVIEW development environment.

LabVIEW includes the following tools to help you analyze your data:

- More than 400 measurement analysis functions for Differential Equations, Optimization, Curve Fitting, Calculus, Linear Algebra, Statistics, etc.
- 12 new Express VIs specifically designed for measurement analysis, including filtering and spectral analysis
- Signal Processing VIs for Filtering, Windowing, Transforms, Peak Detection, Harmonic Analysis, Spectrum Analysis, etc.

Present with LabVIEW




Presentation with LabVIEW can be done on your PC or over a network, or you can take advantage of additional applications such as DIAdem.

LabVIEW includes the following tools to help you present your data:


- On your machine — Graphs, Charts, Tables, Gauges, Meters, Tanks, 3D Controls, Picture Control, 3D Graphs (Windows Only), Report Generation (Windows Only)
- Over the Internet — Web Publishing Tools, Datasocket (Windows Only), TCP/IP, VI Server, Remote Panels, Email
- Enterprise Connectivity Toolset — SQL Tools (Databases), Internet Tools (FTP, Telnet, HTML)

LabVIEW Dialog Box



Project Explorer

- Find, access, and organize project files
- Prevent, detect, and resolve incorrect links
- Deploy or download files to targets
- Manage code for build options
 - Executables, installers, and zip files
- Integrate with source code control providers




LabVIEW Files

- Common LabVIEW file extensions:

LabVIEW project —.lvproj

Virtual instrument (VI) — .vi

Custom control — .ctl




Adding Folders to a Project

– Virtual folder


- Organizes project items and does not represent files on disk

– Auto-populating folder


- Adds a directory on disk to the project
- LabVIEW continuously monitors and updates the folder according to changes made in the project and on disk



Creating a new VI



- Blank VI or File»New VI to open a blank VI




- VI from Template... or File»New... to open the Template Browser ,and configure a VI template, global variable, control, etc...

Parts of a VI – Front Panel

Front Panel – User interface for the VI


You build the front panel with controls (inputs) and indicators (outputs).




Parts of a VI – Block Diagram

Block Diagram – Contains the graphical source code


Front panel objects appear as terminals on the block diagram.



Parts of a VI – Icon/Connector Pane



Icon – Graphical representation of a VI



Connector Pane – Map of the inputs and outputs of a VI



Icons and connector panes are necessary to use a VI as a subVI.

- A subVI is a VI that appears on the block diagram of another VI.
- A subVI is similar to a subroutine or function in a text-based programming language.

Course Exercises Notation

Block Diagram

Front Panel



* Comments *
Do not enter these

Name of object




Random Number (0-1) (Numeric subpalette)
Returns a random number between 0 and 1.


Where it can be found

Description

Front Panel Window



Block Diagram Window



Front Panel and Block Diagram Toolbars



Run button

Continuous Run button

Abort button

Pause/Continue button

Font ring

Alignment ring

Distribution ring

Resize ring

Reorder ring

Context Help Button



Additional Buttons on the
Block Diagram Toolbar

- Execution Highlighting button
- Step Into button
- Step Over button
- Step Out button

 Warning indicator

 Enter button

 Broken Run button


Tools Palette

- LabVIEW automatically selects the tool needed
- Available on the front panel and the block diagram
- A tool is a special operating mode of the mouse cursor
- Use the tools to operate and modify front panel and block diagram objects
- To show the tools palette, select
View » Tools Palette




Controls Palette

- Contains the controls and indicators you use to create the front panel.
- Navigate the subpalettes or use the **Search** button to search the Controls palette.



Functions Palette

- Contains the VIs, functions, and constants you use to create the block diagram.
- Navigate the subpalettes or use the **Search** button to search the Functions palette.




Searching for Controls, VIs, and Functions

Ways to find controls, VIs, and functions:

- Search or navigate the palettes.
 - Controls palette
 - Functions palette
- Search by name of object.
 - Quick Drop dialog box
- Search palettes, *LabVIEW Help*, and `ni.com`.
 - Search text box in toolbar


Searching with Quick Drop

- Lets you quickly find controls, functions, VIs, and other items by name.
- Press the <Ctrl-Space> keys to display the Quick Drop dialog box.



Global Search

Use the Search bar in the top right of the front panel and block diagram windows to search palettes, *LabVIEW Help*, and `ni.com`.



Summary

- Virtual instruments (VIs) have three main parts — the front panel, the block diagram, and the icon and connector pane
- The front panel is the user interface of a LabVIEW program and the block diagram is the executable code
- The block diagram contains the graphical source code composed of nodes, terminals, and wires
- Floating Palettes: Tools Palette, Controls Palette (only when Front Panel Window is active), and Functions Palette (only when Block Diagram Window is active)
- There are help utilities including the Context Help Window and LabVIEW Help